

NEO Shield- Series™

Creating New Dimensions for Neonatal & Child Care

Bubble CPAP System



PRODUCT CATALOGUE



www.neokraftmedical.com




NEOKRAFT™
ISO 13485 Certified Company

Bubble CPAP System



NEOPAP - 505



Range of Respiratory Humidifiers



Servo Humidifier



Non-Servo Humidifier with Display



Non-Servo Humidifier

Infant Resuscitator



RESIPUFF 1050

Infant T-Piece resuscitator

Lightweight easy to operate

No power or battery is required



Resipuff 1050 Infant T-Piece resuscitator is a lightweight stand alone unit suitable for delivery wards, Neonatal care units, Operating rooms and Infant transport resuscitation purpose.

- Infant resuscitator will safely inflate the baby's lungs and provide optimum oxygenation by delivering consistent PIP & PEEP with each breath avoiding the risks associated with under or over inflation of uncontrolled pressures
- The patient T-Piece can connect to Neonatal masks or Endotracheal tubes directly.
- Different sizes of face mask which provides flexible in usage for different babies
- Inspiratory high pressure safety cut-off relief setting
- Peak inspiratory pressure setting (PIP)
- Positive end expiratory pressure setting (PEEP)
- Works on low flow from 5 L/m to 15 L/m
- Medical grade monometer to monitor PIP & PEEP delivered

3-in-1 Weighing Scale



*Unique Digital Weighing Solution
for Infant, Child & Adult*

Infant weigh bed
up to 10 Kgs

Folding child seat
up to 40 Kgs

Adult standing
up to 120 Kgs

Digital Monitoring Device



Front View



Back View

* Battery backup for 6 hours



**Infant Digital Weighing Scale
Model IWS-101**



**3-in-1 Digital Weighing Scale
Model MWS-301**

* Technical specifications may be altered by Neokraft without prior notice, accessories shown in the catalogue are not part of the standard equipment

Manufactured by:



NEOKRAFT MEDICAL PRIVATE LIMITED

18, 1st Cross, 11th Main, Tata Nagar, Sahakara Nagar [P],
BANGALORE - 560 092, INDIA.

Phone : (080) 6500 0249

Mobile : +91 9686 78 5555

Email : info@neokraftmedical.com

Website : www.neokraftmedical.com

Channel Partner