

Bili Flux - Neo 200 Sleek




LED Phototherapy


Treat Neonatal Jaundice with Confidence



Estimated Spectral Irradiance



- Treatment irradiance at Skin level up to $60\mu\text{W}/\text{cm}^2/\text{nm}$ at 30cm with wave length 420 to 500nm.
- Digital timer for total lamp usage and patient exposure.
- Only 30 watts power consumption, up to 80% power saving.
- Up to 30% faster serum bilirubin breakdown compared to conventional Phototherapy.
- Optical design ensures uniform light distribution to the patient surface area.



**LED Phototherapy Stand
Model NEO 200 Sleek**



Radiant Warmer
Model : NEO 510



Radiant Warmer
Model : NEO 300



Resuscitation Unit
Model : NEO 310IRU



Double Surface LED Phototherapy
Model : NEO 220



Infant Trolley
Model : NEO 600



3-in-1 Digital Weighing Scale
Model MWS-301



LED Phototherapy Stand
Model NEO 200



Infant Digital Weighing Scale



Phototherapy Eye Band



Infant Oxygen Hood



Manual Resuscitator

* Technical specifications may be altered by Neokraft without prior notice, accessories shown in the catalogue are not part of the standard equipment

Manufactured by:



NEOKRAFT MEDICAL PRIVATE LIMITED

18, 1st Cross, 11th Main, Tata Nagar, Sahakara Nagar [P],
BANGALORE - 560 092, INDIA.

Phone : (080) 6500 0249

Mobile : +91 9686 78 5555

Email : info@neokraftmedical.com

Website : www.neokraftmedical.com

Channel Partner